

NPS-XY-100A & NPS-XY-100A-UHV 100µm x 100µm Super Invar Stage

Nanopositioning Stage



The NPS-XY-100A and the NPS-XY-100A-UHV NanoMechanisms are piezo scanned flexure guided stage with integrated capacitance position sensors capable of sub- nanometer resolution and reproducibilty. The NPS-XY-100A-UHV has been designed to operate down to 10-9 Torr.

Finite element analysis of the flexure guidance mechanisms has reduced parasitic angular motions to less than 25 micro-radians over the full 100µm range.

The Super Invar construction (CTE 0.3ppm K-1, compared to Aluminium at 23ppm K-1) minimizes thermal drift, which on a nanometer scale can be very important.

The unique iso-static mounting system ensures that stresses from the mounting system are properly relieved and establishes the center of the stage as the co-ordinate reference point.

Key features

- >100µm travel in each axis with sub-nanometer resolution
- Typically <0.005% hysteresis and <0.01% linearity error
- First resonant frequency >350Hz
- In-situ scanning and stepping response optimization
- Robust and reliable
- Super Invar construction

Typical applications

- AFM, SPM, NSOM
- High Precision Microscopy

Suggested controllers

NanoScan NPC-D-6330 Series Controllers Delivering Low noise, Low drift, High power, High Resolution and High position update rates.

Designed specifically to control Queensgate's Nanometer Precision Mechanisms incorporating capacitive sensors. They give precise positional feedback delivering high resolution and linearity of movement. The fast update rate and Queensgate control algorithms contribute to high speed positioning accuracy for dynamic applications that require high speed movement of the stage.

The PC software facilitates user optimisation of all operating parameters, including PID and notch filter set up. There are eight programmable slots, three which are populated to provide fast, medium and slow PID settings, the addition five slots are available for application specific settings.

Function playback provides user defined pre-programmed waveforms for applications such as rasta scanning or constant velocity scanning. The calibration and dynamic settings are held in the stage eprom which allows controllers to be interchanged with minimal performance changes.

NPS-XY-100A & NPS-XY-100A-UHV Nanopositioning Stage

Technical specification NPS-XY-100A

Parameter		Symbol	Value			Units	Comments
Static physical							
Material		Super Invar (Bright nickel plated)					
Size			100 x 100 x 23			mm	Note 1
			Minimum	Typical	Maximum		
*Range closed loop		dxp∙max	± 50	± 55		μm	
*Scale factor error (1σ)		δbx1			0.1	%	
Static stiffness				1		N∙µm-1	
Resonant frequency:	0g load	F0·0		350		Hz	
	50g load	F0·50		260		Hz	
	1000g load	F0·100		120		Hz	
Maximum load					1	Kg	Note 2
Dynamic physical (Typ	oical values)					
			Fast	Medium	Slow		Note 3
*3dB Bandwidth		Bx∙p	53	20	4		Typical
*Small signal settle time		txs·s	15	30	130	ms	Note 4
*Position noise (1σ)		δxp∙n	0.5	0.4	0.25	nmRms	Note 5
Error terms							
			Minumum	Typical	Maximum		
*Hysteresis (peak to peak)		δxp∙hyst		0.005	0.015	%	Note 6
*Linearity error (peak)		δxp·lin		0.005	0.01	%	Note 6
*Rotational error		δφχ		10	25	µradians	Note 7
*Rotational error		δθχ		5	10	µradians	Note 7
*Rotational error		δγχ		5	10	µradians	Note 7
Orthogonality		δθorth		8	10	mradians	

Notes

*These parameters are measured and supplied with each mechanism

- 1. With 40mm diameter central aperture.
- 2. Depends on orientation. 1 kg is the maximum load for gravity acting in the Z direction. 0.5 kg is the maximum load for gravity acting in the X or Y axes. Loads greater than 5 kg can cause damage to the flexure mechanism.
- 3. For dynamic operation the servo-loop parameters are preset for different performances; the parameters are user settable via software control. Fast means the fastest the stage can stably move with less than 50 grams load. Medium means the maximum stable speed for loads up to 200 grams. Slow means the speed at which the servo loop is stable for all masses up to the maximum allowed mass equivalent to low noise setting.
- 4. Step and Settle time is the time taken to settle to within 2% of the step measured using an interferometer. The step settle time is a function of the servo loop parameters which are user controllable. The test step size is 2000 nm.
- 5. The actual position noise of the stage. Measured with a laser interferometer sampling 1 Hz to 25 kHz.
- 6. Percent of the displacement.
- 7. Percent error over the full range of motion

General Note: Specifications assume standard 2m cable length, for custom cable lengths and UHV position noise, linearity and hysteresis may vary.



WORLDWIDE DISTRIBUTION

Jena, Germany

t: +49 (0) 3641 675 650

e: jena@prior.com



Prior Scientific Ltd Cambridge, UK

t: +44 (0)1223 881711 **e:** uksales@prior.com



Rockland MA U.S.A t: +1 781-878-8442 e: info@prior.com

.	Scientific Gr	





e: info-japan@prior.com



t: +86 (0)512 6617 5866 e: info-china@prior.com



For more information visi